

# **Research Standing Frame**

## **User Guide**



## Research Standing Frame User Guide

### ***Labview Program Location:***



<\\cscharkemafs01\Programs\Research Standing Frame V1.0\ResearchStandingFrame-V1.0.vi>

### ***Labview Version***


Labview 2012

### ***Force Direction***

Top View



Front View



## Research Standing Frame User Guide

### ***BNC Connection***

#### **Input**

Left force plate – force plate 2

Force Plate 2-1 (X12)	→	AI 0	Force Plate 2-2 (X34)	→	AI 1
Force Plate 2-3 (Y12)	→	AI 2	Force Plate 2-4 (Y34)	→	AI 3
Force Plate 2-5 (Z1)	→	AI 4	Force Plate 2-6 (Z2)	→	AI 5
Force Plate 2-7 (Z3)	→	AI 6	Force Plate 2-8 (Z4)	→	AI 7

Right force plate – force plate 3 (BNC cable with green electrical tape)

Force Plate 3-1 (X12)	→	AI 16	Force Plate 3-2 (X34)	→	AI 17
Force Plate 3-3 (Y12)	→	AI 18	Force Plate 3-4 (Y34)	→	AI 19
Force Plate 3-5 (Z1)	→	AI 20	Force Plate 3-6 (Z2)	→	AI 21
Force Plate 3-7 (Z3)	→	AI 22	Force Plate 3-8 (Z4)	→	AI 23

Left load cell (LLC, BNC cable with blue electrical tape)

LLC X → AI 32      LLC Y → AI 33      LLC Z → AI 34

Right load cell (RLC, BNC cable with red electrical tape)



RLC X → AI 35      RLC Y → AI 36      RLC Z → AI 37

DAQ Board Connections	
LX	32
LY	33
LZ	34
RX	35
RY	36
RZ	37
LFX12	0
LFX34	1
LFY14	2
LFY23	3
LFZ1	4
LFZ2	5
LFZ3	6
LFZ4	7
RFX12	16
RFX34	17
RFY14	18
RFY23	19
RFZ1	20
RFZ2	21
RFZ3	22
RFZ4	23

#### **Output**

Counter I/O Ctr0 → Sync pulse output


***Power Switches of Load Cells and Force Plates***




# Research Standing Frame User Guide

## Software Interface

### Force Graphs




### Force VS Time & Baseline





# Research Standing Frame User Guide

## Start the Program

1. Click run button 
2. Remove baseline force reading from load cells and force plates


Click “Baseline” button in “Baseline” tab0




3. Go back to “Force Graph” tab and click “Save” button to save data



# Research Standing Frame User Guide




Two files will be saved (one is force data and one is voltage data).



4. Click “Stop” button to stop program

# Research Standing Frame User Guide



In the dialog box, click “Yes” to stop program and “Cancel” to resume program



### *Data Format*

## Force data

StandingFrameNewtonTestData - Notepad

File Edit Format View Help


STANDING FRAME NEWTON FORCE DATA  
SCAN\_RATE: 1000 HZ  
SAMPLE\_PER\_PULSE: 1

LX, LY, LZ, RX, RY, RZ, LF X12, LF X34, LF Y14, LF Y23, LF Z1, LF Z2, LF Z3, LF Z4, RF X12, RF X34, RF Y14, RF Y23, RF Z1, RF Z2, RF Z3, RF Z4 (UNIT: Newton)

0.002, 0.027, 0.012, 0.013, 0.002, 0.002, -0.006, -0.006, -0.022, -0.126, -0.157, 0.008, 0.026, 0.003, 0.147, 0.015, 0.073, -0.063, 0.085, -0.000, -0.044, 0.152  
0.004, -0.008, 0.017, 0.010, 0.021, 0.018, -0.109, 0.020, 0.004, -0.074, -0.183, -0.043, 0.051, 0.003, 0.017, 0.002, 0.060, -0.167, -0.119, -0.000, -0.044, 0.076  
0.004, 0.013, 0.007, -0.014, 0.016, 0.002, 0.046, 0.020, 0.069, -0.087, 0.020, 0.033, 0.076, 0.003, 0.082, -0.011, 0.124, -0.024, 0.008, -0.000, -0.044, 0.127  
0.007, 0.019, -0.010, -0.014, -0.001, 0.024, -0.109, 0.033, -0.022, -0.178, -0.030, -0.018, 0.026, 0.003, 0.069, -0.128, -0.005, 0.001, -0.042, 0.025, 0.006, 0.203  
0.002, 0.008, 0.001, -0.001, -0.001, -0.009, 0.072, -0.110, 0.004, -0.074, -0.056, 0.033, 0.076, 0.003, 0.043, -0.037, 0.060, -0.128, -0.017, -0.000, -0.044, -0.101  
0.002, -0.008, 0.017, -0.003, 0.016, 0.029, 0.124, 0.033, -0.073, -0.139, 0.097, 0.059, 0.026, 0.003, 0.160, 0.028, 0.047, -0.063, -0.220, -0.000, 0.006, 0.076  
0.007, 0.016, 0.001, 0.005, 0.007, 0.035, -0.187, 0.007, -0.061, -0.190, -0.107, 0.008, 0.051, -0.022, 0.004, 0.041, 0.047, -0.011, 0.034, -0.000, -0.019, -0.202  
-0.001, 0.003, 0.017, -0.009, 0.002, 0.013, -0.031, -0.058, 0.009, -0.126, -0.284, 0.033, 0.051, -0.022, 0.056, -0.037, 0.124, -0.102, -0.017, -0.000, 0.006, -0.126  
0.004, 0.024, 0.001, -0.003, 0.013, -0.009, -0.264, -0.045, 0.017, -0.307, 0.020, -0.018, -0.049, 0.003, 0.082, 0.054, 0.111, -0.063, -0.042, -0.000, -0.019, 0.178

Ln1, Col1

## Voltage data



StandingFrameVoltageTest6.dat - Notepad

File Edit Format View Help

STANDING FRAME VOLTAGE DATA

SCALING RATE: 1000 HZ

SAMPLE PER PULSE: 1

LX, LY, LZ, RX, RY, RZ, LF L2 (UNIT: voltage)

0.073,	0.019,	0.137,	-0.007,	0.049,	0.117,	0.004,	0.004,	0.002,	0.009,	0.009,	0.007,	0.005,	0.010,	0.004,	0.014,	0.002,	0.009,	0.006,	0.010,	0.010,	0.005
0.073,	0.019,	0.137,	-0.007,	0.049,	0.117,	0.004,	0.004,	0.002,	0.009,	0.009,	0.007,	0.005,	0.010,	0.003,	0.014,	0.002,	0.009,	0.005,	0.010,	0.010,	0.005
0.073,	0.019,	0.137,	-0.007,	0.049,	0.117,	0.004,	0.004,	0.002,	0.009,	0.009,	0.007,	0.005,	0.010,	0.004,	0.014,	0.002,	0.009,	0.006,	0.010,	0.010,	0.005
0.073,	0.019,	0.137,	-0.007,	0.049,	0.117,	0.004,	0.004,	0.002,	0.009,	0.009,	0.007,	0.005,	0.010,	0.004,	0.014,	0.002,	0.009,	0.005,	0.010,	0.010,	0.005
0.073,	0.019,	0.137,	-0.007,	0.049,	0.117,	0.004,	0.004,	0.002,	0.009,	0.009,	0.007,	0.005,	0.010,	0.004,	0.014,	0.002,	0.009,	0.005,	0.010,	0.010,	0.005
0.073,	0.019,	0.137,	-0.007,	0.049,	0.117,	0.004,	0.004,	0.002,	0.009,	0.009,	0.007,	0.005,	0.010,	0.004,	0.014,	0.002,	0.009,	0.005,	0.010,	0.010,	0.005
0.073,	0.019,	0.137,	-0.007,	0.049,	0.117,	0.004,	0.004,	0.002,	0.009,	0.009,	0.007,	0.005,	0.010,	0.004,	0.014,	0.002,	0.009,	0.005,	0.010,	0.010,	0.005
0.073,	0.019,	0.137,	-0.007,	0.049,	0.117,	0.004,	0.004,	0.002,	0.009,	0.009,	0.007,	0.005,	0.010,	0.004,	0.014,	0.002,	0.009,	0.005,	0.010,	0.010,	0.005
0.073,	0.019,	0.137,	-0.007,	0.049,	0.117,	0.004,	0.004,	0.002,	0.009,	0.009,	0.007,	0.005,	0.010,	0.004,	0.014,	0.002,	0.009,	0.005,	0.010,	0.010,	0.005
0.073,	0.019,	0.137,	-0.007,	0.049,	0.117,	0.004,	0.004,	0.002,	0.009,	0.009,	0.007,	0.005,	0.010,	0.004,	0.014,	0.002,	0.009,	0.005,	0.010,	0.010,	0.005
0.073,	0.019,	0.137,	-0.007,	0.049,	0.117,	0.004,	0.004,	0.002,	0.009,	0.009,	0.007,	0.005,	0.010,	0.004,	0.014,	0.002,	0.009,	0.005,	0.010,	0.010,	0.005

BASELINE VALUES

LX, LY, LZ, RX, RY, RZ, LF L2

0.073,	0.019,	0.137,	-0.007,	0.049,	0.117,	0.004,	0.004,	0.002,	0.009,	0.009,	0.007,	0.005,	0.010,	0.003,	0.014,	0.002,	0.009,	0.006,	0.010,	0.010,	0.005
--------	--------	--------	---------	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	-------

Ln1, Col1